

## ID Screen® PPR Competition

Competitive ELISA for the detection of anti-PPRV nucleoprotein antibodies

Peste des petits ruminants (PPR) is a contagious disease affecting goats and sheep primarily in Africa, the Middle-East and the Indian subcontinent. It is caused by a species of the *Morbillivirus* genus of viruses. The disease is highly contagious, with approximately 80 percent mortality in acute cases. In June 2008, the disease invaded Morocco, crossing the natural barrier of Sahara and causing concern that the disease could spread into Europe.

Serology may be used to identify and control outbreaks. The ID Screen® PPR Competition ELISA efficiently detects antibodies directed against the virus nucleoprotein.

All components are ready-to-use and each sample is deposited only once.

The test uses technology developed by the OIE reference laboratory (CIRAD-EMVT, Montpellier, France).

Product Code	PPRC-4P
Test Principle	Competitive ELISA
Antigen	PPR recombinant nucleoprotein
Conjugate	Anti-NP-HRP concentrated conjugate (10)
Sample Type	Multiple species, including sheep and goat serum and plasma. (Please contact IDVET for more information.)
Sample dilution factor	1:2
Number of tests	384 (4 plates); Please contact IDVET for other formats
Microplate format	12 x 8-well strips
Protocol	<ol style="list-style-type: none"> <li>1. Sample Incubation 45 min</li> <li>2. Conjugate Incubation 30 min</li> <li>3. Three washes</li> <li>4. Substrate Incubation 15 min</li> </ol>
Test Interpretation	<p>S/N &lt; 50% = positive</p> <p>50 % &lt; S/N ≤ 60% = doubtful</p> <p>S/N &gt; 60% = negative</p>

### Reference

Development of a competitive ELISA for detecting antibodies to the Peste des Petits Ruminants virus using a recombinant nucleoprotein. Libeau G, Préhaud C, Lancelot R, Colas F, Guerre L, Bishop DH, Diallo A., Res Vet Sci. 1995 Jan;58(1):50-5.




Screening format:  
each sample is  
deposited only once

Ready-to-use  
components,  
including coated  
plates

Simple and easy-to-  
use: results in 90  
minutes

OIE reference lab  
technique

High specificity and  
sensitivity



### Kit Contents

Coated microplates  
Concentrated Conjugate (10X)  
Positive Control  
Negative Control  
Dilution Buffers  
Wash Concentrate (20X)  
Substrate Solution (TMB)  
Stop Solution